

# ambIT® PIB•PCA

## Superior Pain Management

### PIB•PCA Infusion Pump

Better Distribution of Medication  
More Effective Pain Management

#### Three Unique Programs:

- **PCA** (Patient Controlled Analgesia) with continuous rate and optional patient-controlled bolus
- **PIB** (Programmable Intermittent Bolus) with regularly occurring programmed doses
- **PIB + PCA** (Combined PIB and PCA) with regularly occurring doses and manual bolus options

**ambIT® PIB•PCA is designed to meet patients' needs by combining qualities of PCA and intermittent systems into a single pump**

**ADMEDUS**



#### Clinical Applications

- Continuous Regional Nerve Blocks
- Intermittent Bolus Nerve Blocks
- Surgical Site Infusions

# ambIT® PIB•PCA

Volumetric Accuracy	± 6 %
Basal Flow Rates	0.1 to 20 ml/hr in 0.1 ml/hr steps
Bolus Delivery Rate	210 ml/hr
Bolus Dose Volume	0 to 50 ml in 0.1ml steps
Bolus Lockout Times	00:01 to 24:00 in 1 minute increments
Volumes to be Infused (ml)	0 to 1000
Alarms & Alerts	<ul style="list-style-type: none"> <li>• Detached Cassette</li> <li>• Malfunction</li> <li>• Run Indicator Light</li> <li>• Bolus Infusing</li> <li>• Downstream Occlusion</li> </ul>
	<ul style="list-style-type: none"> <li>• Low Battery</li> <li>• Dead Battery</li> <li>• Bolus Requested</li> <li>• Infusion Complete</li> <li>• Pause Indicator</li> </ul>
	<ul style="list-style-type: none"> <li>• Volume Infused</li> <li>• Boluses Requested</li> <li>• Doses Requested</li> </ul>
	<ul style="list-style-type: none"> <li>• Boluses Delivered</li> <li>• Time Elapsed</li> </ul>
	<ul style="list-style-type: none"> <li>• First dose can be delivered immediately or after delay</li> </ul>
Reports	<ul style="list-style-type: none"> <li>• Volume Infused</li> <li>• Boluses Requested</li> <li>• Doses Requested</li> </ul>



## ambIT® PIB•PCA Infusion Pump and Accessories

Product description	Code
ambIT® PIB•PCA Infusion Reusable Pump	220501
ambIT® PIB•PCA Single Use Kit - Pump, Luer Lock Cassette, 420ml Medibag, Patient Carry Bag, Button Guard, Peel Away Introducer 17 Gauge 203mm, 2 'AA' Batteries. Catheter with 127mm fenestrations	220507
ambIT® PIB•PCA Single Use Kit - as above. Catheter with 254mm fenestrations	220508
ambIT® PIB•PCA Single Use Kit - as above. Catheter with 76mm fenestrations	220509
ambIT® PIB•PCA Single Use Kit Pump, Spike Cassette, Patient Carry Bag, Button Guard, 2 'AA' Batteries	220513
ambIT® Cassette, Filter, Spike - also available with yellow stripe (Y)	220140 220514-Y
ambIT® Cassette, luer lock for use with Medibag - also available with yellow stripe (Y)	220266 220266-Y
ambIT® Medibag 420ml	220278