

# ambIT<sup>®</sup> PreSet PCRA

When effective pain control is paramount

## Patient Controlled Regional Anaesthesia

Accurate Safe Mobile Easy to use

- Volumetric accuracy of +/- 6%
- Choice of three delivery profiles
  - basal rate and bolus
  - basal rate only
  - bolus only
- Bolus lockout period
- Program can be adjusted during treatment
- Four preset and one individualised program
- History reports available



### Clinical Applications:

- Continuous Regional Nerve Blocks
- Surgical Site Infusions
- Epidural Infusions
- IV PCA
- Sub-cutaneous Infusions

**ADMEDUS**

#### Customer Service

Australia T 1300 550 310 E info.au@admedus.com  
New Zealand T 0800 226 926 E info.nz@admedus.com

admedus.com

# It's all about accuracy.

## Product Comparison

### Elastomeric Pumps

- Limited programming flexibility
- Less accurate flow-rate
- No alerts or alarms
- High-cost for low technology



### ambIT® PreSet PCRA

- Simple, programmable operation
- Volumetric accuracy of  $\pm 6\%$
- Visual and audible safety alerts
- Low-cost for high technology



Avoid the variability of elastomeric technology

## ambIT® PreSet™ PCRA single use kit optional components:

- ambIT® pump
- Luer lock cassette
- Spike cassette
- 420ml medi-bag
- 2 AA batteries
- Patient carry bag
- Button guard
- Insertion trocar with peel away sheath
- Surgical site catheter (choice of 76mm, 127mm or 254 fenestrations)



Please contact us to discuss your kit requirements.