


# DOSI-FUSER<sup>®</sup>

## Advanced Elastomeric Infuser

### PORTABLE ELASTOMERIC PUMP

#### FEATURES

- **Ambulatory** - allows patients to fully mobilise
- **No complex programming**
- **Safe and accurate** - internal gauge
- **Full range of flow rates**



#### Clinical Applications

- Acute and chronic pain management
- Surgical site and nerve block infusions
- Antibiotic and chemotherapy delivery
- Hydration

# ADMEDUS

#### Customer Service

Australia T 1300 550 310  
New Zealand T 0800 226 926

[admedus.com](http://admedus.com)

## DOSI-PAIN® KIT

Includes:

- Elastomeric pump with one or two channels
- One or two DOSI-PAIN® sets (catheter, introducer, syringe and dressings)
- 50ml Luer-lock syringe for filling pump
- Mesh carry bag



Catheter fenestrations 25mm, 80mm, 160mm, 240mm and 300mm.

## DOSI-FUSER® High Volume

Configuration:

- 400ml from > 1 hour - 2 days
- 500ml from > 1 hour - 2 days
- 600ml from > 1 hour - 2 days



## DOSI-FUSER® Multi-Flow

- Seven gradually increasing flow rates
- Allows physician to customise patients pain protocol
- Safe lock with removable key

Configurations:

- 100ml: 1-7ml/hr, 2-14 ml/hr
- 250ml: 1-7ml/hr, 2-14 ml/hr

